

Protocol 0370

## ***ElectroSquarePorator™***

T820/830 Electroporation Protocol

**Cell line: HepG2, Human Hepatocellular Carcinoma (ATCC HB 8065)**

Transfectant: pGL3 (Promega SV40 promoter & enhancer, luciferase)

---

### **Cell Preparation:**

Growth medium: Grow cells in DME/F12 containing antibiotics and 10% FBS. Cells must be in good condition (exponential growth).

Trypsinization: Wash with FBS free media, add trypsin; inactivate trypsin with 10% FBS.

Washing Procedure: Wash cells 2x using DME/F12 media with 2.5% FBS (no antibiotics)

Cell Density: Resuspend @  $1 \times 10^7$  cells/ml.

---

### **Electroporation Settings:**

Choose Mode: LV Mode (99 msec/500V)

Set Voltage: 150 V

Set Pulse Length: 70 msec

Set Number of Pulses: 1

Chamber: BTX Disposable Cuvette Model 640, 4mm gap

Desired Field Strength: 375 V/cm

---

### **Electroporation Procedure:**

Sample Volume: 400  $\mu$ l

DNA Concentration: 50  $\mu$ g/ml, or 20  $\mu$ g total

Temperature: Room temperature (5 min. pre and post incubation)

Pulse: Press the Start button to activate the Automatic Charge & Pulse sequence

Post Pulse Treatment: Plate  $0.4 \times 10^6$  transfected cells into a single well of a 24 well plate. Add 1 ml of DME/F12 media with 10% FBS and incubate for 24 hours @ 37°C.

---

**Results:** 80-90 % cell viability. Luciferase activity 1,700x above background.

**Reference:** Personal communication, Dr. David Pasco, National Center for the Development of Natural Products, University of Mississippi