

**Protocol 0906****ElectroSquarePorator<sup>™</sup>****ECM<sup>®</sup> 630 HT 25/96 Multi-Well ELECTROPORATION PROTOCOL**Cell Line: *Sacchomyces cerevisiae*

Transfectant: plasmid DNA

**Cell Preparation:** *Sacchomyces cerevisiae* were grown overnight to stationary phase (OD<sub>600</sub> of ~3). An aliquot of the culture was inoculated into 100 ml of YPD media (10 g/l yeast nitrogen base, 20 g/l peptone and 20 g/l D-(+)-glucose) to reach OD<sub>600</sub> of ~0.3. The cells were grown until OD<sub>600</sub> reached ~1.6 before collecting by centrifugation. The cell pellet was washed twice with 50 ml of cold water, and once with 50 ml electroporation buffer (1 M sorbitol/1 mM CaCl<sub>2</sub>). Cells were conditioned in 20 ml of 0.1 M LiAc/10mM DTT by incubation in culture flask with shaking (225 rpm) for 30 min at 30°C, washed one or more times with 50 ml of electroporation buffer, then suspended in 100-200 µl of the same buffer to reach 1 ml volume, which corresponds to 1.6 x 10<sup>9</sup> cells/ml. (The scale up of culture and washed volumes should be proportional if more electrocompetent yeast suspension is needed.)

**Electroporation Settings:**


Choose Mode:	HV
Set Voltage:	2500 V
Electrode Gap:	2 mm*
Capacitance:	25 µF
Resistance:	~100 – 200 Ω
Pulse number:	1
Electrode Type:	96 – well plate
Desired Field Strength:	12.5 kV/cm
Time constant:	3 – 4.5 ms

**Electroporation Procedure:**

Temperature:	4°C
Total Volume:	250 µl
Transfectant conc:	4 µg vector:12 µg insert (1:3 ratio)

**Many Challenges,  
One Solution.**[www.btxonline.com](http://www.btxonline.com)

Cell density:  $1.6 \times 10^9$  cells/mL  
Pulse: Press Automatic **Start** to activate Charge and Pulse Sequence  
Post Pulse Treatment: Cells were suspended in 10 ml of 1:1 mix of 1 M sorbitol: YPD media and incubated at 30°C for 1 h.



## Results:

Increased DNA input and high voltage, but not vector to insert ratio, are critical for maximal transformation efficiency. *Saccharomyces cerevisiae* were treated as described in Fig. 1. Four hundred microliters of yeast cell suspension ( $1.6 \times 10^9$  cells/mL) were electroporated using 1, 2, 3 or 4 mg vector DNA per cuvette (vector to insert ratio is always maintained at 1:3). After 72 h the numbers of colonies were determined and represented as (A) transformation efficiency, transformants/mg vector DNA (mean+standard deviation) and (B) transformation yield, total number of transformants per cuvette (mean+ standard deviation). Additionally, (C) 400 µl yeast suspension were electroporated by 2 or 2.5 kV and different vector to insert DNA ratios. The transformation efficiencies are expressed as in (A).

**Reference:** Benatui, L. et al., An improved yeast transformation method for the generation of very large human antibody libraries, February 3, 2010, Protein Engineering, Design & Selection vol. 23 no. 4 pp. 155-159

### BTX item numbers:

450423 ECM 630 generator, 2 mm gap 2 x 96-well plates, Plate seals and HT-200 plate handler

\*Protocol adapted for a 96-well plate.

**Many Challenges,  
One Solution.**

[www.btxonline.com](http://www.btxonline.com)

Compiled: November 1, 2010 – DMerdaa, BTX Technical Support