

## Protocol 1040

# GEMINI & ECM® 830 ELECTROPORATION PROTOCOL

**Cell Line:** Primary Human CD4 and CD8 T lymphocytes (CAR T cells)

**Application or Transfectant:** CRISPR - Cas9 mRNA

**Electrode:** 2 mm gap cuvette Item # 45-0125

---

### Cell Preparation:

Primary human CD4 and CD8 T cells were isolated from healthy volunteer donors following leukapheresis. The T cells were washed three times with OPTI-MEM and re-suspended in OPTI-MEM (Invitrogen) at a final concentration of  $1$  to  $3 \times 10^8$  cells/ml. The recommended electroporation buffer with the current ECM 830 model is BTXpress high performance buffer. Electroporation of CRISPR reagents with one-shot CAR or CAR T cells was performed with a BTX ECM 830 electroporator.

---

### Square Wave Electroporation Settings:

|                         |           |
|-------------------------|-----------|
| Set Voltage:            | 360 V     |
| Set Pulse Length:       | 1 ms      |
| Set Number of Pulses:   | 1         |
| Desired Field Strength: | 1800 V/cm |

---

### Electroporation Procedure:

|                         |   |
|-------------------------|---|
| Electrode Gap:          | 2 mm gap cuvette, Item # 45-0125  |
| Total Sample Volume:    | 100 $\mu$ l   |
| Number of Cells:        | $1$ to $3 \times 10^8$ cells/ml $\times 10^6$ PBMCs/ml  |
| Electroporation Buffer: | Opti-MEM or BTXpress High Performance Electroporation Solution (Item # 45-0802)   |
| Amount of Transfectant: | 20 $\mu$ g cas9 RNA (per TCR chain)   |
| Pulse:                  | Press the Go icon or the Start button to activate the automatic charge and pulse sequence   |
| Post Treatment:         | The cells were immediately placed in 2 ml of pre-warmed culture media and cultured in the presence of IL-2 (100 IU/ml) at 37°C and 5% CO <sub>2</sub> . |

---

### Results:

The research group developed an easy-to-use immunoassay kit to control T cell bioassays by transfecting RNA-encoded TERS (TCR-engineered reference samples) into primary lymphocytes. High transfection efficiencies of up to 97.3% antigen-specific CD8+ T cells and high viability up to 96.0% were obtained with the BTX ECM 830 Electroporation System.

**Efficiency:** 97.3% antigen-specific CD8+ T cells

**Viability:** 96%




Figure 2 from Bidmon et al., Optimized electroporation settings for six commonly used electroporation devices including: the BTX ECM 830 from Harvard Bioscience; the Eppendorf – Multiporator; the Bio-Rad - Genepulser II and Xcell; and, the Lonza - Nucleofector II and Amaxa 4D. The BTX ECM 830 achieved TCR expression levels higher than the benchmark of > 80% TCR+ of CD8+ and a viability higher than the benchmark > 90% viability.

## References:

Bidmon N et al., [Development of an RNA-based kit for easy generation of TCR-engineered lymphocytes to control T-cell assay performance](#). J. Immunol. Methods 2018; 458: 74–82.

Compiled by: Radhika Dixit, PhD, BTX Applications Scientist, 10/17/2018