

# Tech-Trends

## Volume 3, Series 4

### siRNA transfection of human mesenchymal progenitor cells by electroporation

#### Activin A expression regulates multipotency of mesenchymal progenitor cells

##### Background:

Bone marrow (BM) stroma currently represents the most common and investigated source of mesenchymal progenitor cells (MPCs); however, comparable adult progenitor cells or stem cells have also been isolated from a wide variety of tissues. This study aims to assess the functional similarities of MPCs from different tissues and to identify factor(s) related to their multipotency.

##### Methods:

For this purpose, we directly compared MPCs isolated from different adult tissues, including bone marrow, tonsil, muscle, and dental pulp. We first examined and compared proliferation rates, immunomodulatory properties, and multidifferentiation potential of these MPCs *in vitro*. Next, we specifically evaluated activin A expression profile and activin A: follistatin ratio in MPCs from four sources.

##### Results:

The multidifferentiation potential of the MPCs is correlated with activin A level and/or the activin A: follistatin ratio. Interestingly, by siRNA-mediated activin A knockdown, activin A was shown to be required for the chondrogenic and osteogenic differentiation of MPCs. These findings strongly suggest that activin A has a pivotal differentiation-related role in the early stages of chondrogenesis and osteogenesis while inhibiting adipogenesis of MPCs.

##### Conclusions:

This comparative analysis of MPCs from different tissue sources also identifies bone marrow-derived MPCs as the most potent MPCs in terms of multilineage differentiation and immunosuppression, two key requirements in cell-based regenerative medicine. In addition, this study implicates the significance of activin A as a functional marker of MPC identity.




**ECM 830 Square Wave Electroporation System includes ECM 830 Generator, 630B Safety Stand, Cuvettes 1 mm, 2 mm and 4 mm pkg. of 30 and Cuvette Rack 660**

Cat. 450002

##### References:

Fatès Drouot<sup>1,2</sup>, Wesley M. Jackson<sup>1</sup>, Brent Bobick<sup>1,3</sup>, Sasa Janjani<sup>1,4</sup>, Yingjie Song<sup>1</sup>, George T.J. Huang<sup>1,5</sup> and Roddy S. Tuan<sup>1,6</sup>  
<sup>1</sup>Cartilage Biology and Orthopaedics Branch, National Institute of Arthritis and Musculoskeletal and Skin Diseases, National Institutes of Health, Department of Health and Human Services, 9000 Rockville Pike, Bethesda, Maryland 20892, USA; <sup>2</sup>Current address: Inserm U 844, Hôpital Saint-Eloi, Bâtiment RM, 34295 Montpellier, France; <sup>3</sup>Current address: Department of Biological Sciences, 2500 University Drive NW, Calgary, Alberta, T2N 1N4, Canada; <sup>4</sup>Current address: Department of Otorhinolaryngology, Head & Neck Surgery, Zagreb Clinical Hospital Center, Zagreb University School of Medicine, Salata 4, 10000 Zagreb, Croatia; <sup>5</sup>Current address: Department of Molecular and Cell Biology, Boston University Henry M. Goldman School of Dental Medicine, 72 East Concord Street, Boston, MA 02118, USA; <sup>6</sup>Current address: Department of Orthopaedic Surgery, Center for Cellular and Molecular Engineering, University of Pittsburgh School of Medicine, 450 Technology Drive, Room 221, Pittsburgh, PA 15219, USA

Stem Cell Research & Therapy 2010, 1:11  
<http://stemcellres.com/content/1/2/11>



**Role of activin A in MPC multipotency.** BM-MPCs transiently transfected either with control siRNA or with activin A siRNA were analyzed. (a) Activin A secretion level in the supernatant of transfected BM-MPCs quantified by ELISA 3 days (D3) after transfection showed the efficiency of knockdown by using activin A siRNA. (b, c) Effect of activin A knockdown on BM-MPC and T-MPC chondrogenic potential. Chondrogenic differentiation was evaluated after 21 days in micro pellet culture. (b) Immunohistochemical analysis of cartilage markers (COL2 and AGN) in pellet cultures of BM-MPCs and T-MPCs as well as alizarin blue (AB) and vonkries red (PS) staining revealed that, compared with the untransfected and Control siRNA, MPCs transfected with activin A siRNA exhibited less chondrogenic potential. (c) Quantitative RT-PCR also showed decreased expression of COL2 in MPCs transfected with activin A siRNA, compared with the cells transfected with control siRNA (c). (d, e) Effect of activin A knockdown on BM-MPC osteogenic potential. Both the expression levels of the osteogenic markers (OC and ALP), determined by quantitative RT-PCR, and alizarin red staining, revealed a decrease in the osteogenic activity of BM-MPCs with activin A knockdown. (b, c) Effect of activin A knockdown on BM-MPC adipogenic potential. Both expression levels of the adipogenic markers (PPAR-γ and LPL), determined by quantitative RT-PCR and of oil red O staining, revealed an increase in the adipogenic potential of BM-MPCs with activin A knockdown, compared with BM-MPCs transfected with the irrelevant Control siRNA. RT-PCR data represent the mean ± SD of three independent experiments. \*P < 0.05, versus Control siRNA.

#### ECM 830 ELECTROPORATION PROTOCOL

##### Cell Preparation:

For adipogenic differentiation, cells were seeded at a density of 20,000 cells/cm<sup>2</sup> and treated for 3 weeks with adipogenic medium consisting of DMEM with 10% FBS, and supplemented with 0.5 mmol/L 3-isobutyl-1-methyl-xanthine (IBMX), 1 μg/ml insulin, and 1 μmol/L dexamethasone. For osteogenic differentiation, cells were seeded into 6-well plates at a density of 20,000 cells/cm<sup>2</sup>, and treated for 3 weeks with osteogenic medium, consisting of DMEM with 10% FBS, and supplemented with 10 mmol/L β-glycerol-phosphate, 10 nmol/L dexamethasone, 50 μg/ml ascorbic acid-2-phosphate, and 10 nmol/L 1,25 di-hydroxyvitamin D<sub>3</sub>. To induce chondrogenic differentiation, 96-microwell polypropylene plates were seeded with a density of 300,000 cells per well, and cell pellets formed by centrifugation at 1100 rpm for 6 minutes. The pellet cultures were treated for 3 weeks with chondrogenic medium, consisting of high-glucose DMEM supplemented with 100 nmol/L dexamethasone, 50 μg/ml ascorbic acid-2-phosphate, 100 μg/ml sodium pyruvate, 40 μg/ml L-proline, 10 ng/ml recombinant human transforming growth factor-β<sub>3</sub>, and 50 mg/ml ITS-premix stock.

##### Electroporation settings:

Voltage: 1075 V  
 Pulse Length: 105 μs  
 Pulse Interval: 5 s  
 Pulse Number: 2  
 Electrode Gap: 2 mm  
 Pulse: Press Start to Activate automatic charge and pulse  
 Post treatment: Incubate for 15 min at 37°C in a 5% CO<sub>2</sub>/95% air atmosphere.

##### Electroporation procedure:

Total volume: 300 ul  
 Cell density: 750,000  
 Transfectant: 45 ul of 22.5 ug siRNA



Molecular Delivery Systems